

CARE4 MODULE NAME: WEB Acute/ED/AMB

My Widget Customization

Overview

This document provides instructions on how to view and customize My Widget.

My Widget can be customised to include patient data, such as patient's weight, specific vitals, specific LABs, and reports for a patient's visit.

My Widget can be accessed in the Reference Region of the patient's Chart section.

The screenshot shows the CARE4 patient chart interface for a patient named Sally Test. The interface is divided into several sections: Diagnostics, Provider Notes, Nurse/Allied Health, and Medications. The My Widget section is highlighted with a red box, indicating it is the focus of the document. The My Widget section currently displays "No Data to Display".

Adding Items to My Widget

- Pulling information from other widgets

Right click on the header for the item you want to add to My Widget

Ensure **My** is bolded in the overlay


Click on **My Widget** label to add content to My Widget

NOTE: If **Patient** is bolded it will add the content to the Patient Widget and not My Widget


The screenshot shows the context menu for the "Blood Type (Auto)" widget. The menu options are: Graph, Graph Favorites (Blood Pressure, Pulse Rate, Temperature, Weight, BMI), My Patient Widget (highlighted with a red box and a red arrow), and Queue for Reference. A red box also highlights the "Blood Type (Auto)" header in the background. A red arrow points to the "Click Here" text in the menu.

➤ Pulling information using the Search field in the Reference Region

Enter the query name you would like to pull information from - e.g., weight
Click on the query listed in the drop down list



On the overlay Right Click the result
Ensure **My** is bolded in the overlay
Click on **My Widget** label to add content to My Widget



- Pulling information from other tabs within the patients EMR

E.g., pulling WBC LAB into My Widget

Click on Diagnostics EMR tab
 Click Laboratory tab
 Right Click on WBC

The screenshot shows the EMR interface with the 'Diagnostics' tab selected. Below it, the 'Laboratory' sub-tab is highlighted. A table of lab results is displayed with the following data:

	05/09/2021 16:03	01/10/2021 13:29	22/10/2021 14:12
WBC (4.0-11.0 x10 ⁹ /L)	5.0	7.71	
RBC (4.50-6.50 x10 ¹² /L)	4.70	3.68 L	
Hgb (135-180 g/L)	120 L	109 L	

- Data will now be available through My Widget

The 'My Widget' section displays the following data points:

Weight	0 g	13/11/2021
WBC (4.0-11.0)	7.71 x10 ⁹ /L	01/10/2021

Removing/Editing My Widget Items

Right Click item within the My Widget
On the overlay click **EDIT** for **My Widget**

The screenshot shows a right-click context menu for a widget. The menu items are: Graph, Graph Favorites (with sub-items: Blood Pressure, Pulse Rate, Temperature, Weight, BMI), My | Patient Widget (with sub-items: Weight, White Blood Count), Add to My Flowsheet, and Queue for Reference. The 'Edit' button next to 'My | Patient Widget' is circled in red. The widget itself is titled 'External Application' and 'My Widget'. It displays a table with one row: 'Weight' 120 kg 07/09/2021. The 'Weight' column is also circled in red. Below the table is a 'Patient Widget' section with 'No Data to Display'.

The My Widget overlay will display
Sort queries using the up and down arrows ↓ ↓
Remove queries by clicking the **X**

This screenshot shows the 'My Widget' overlay menu with sorting arrows (up and down) and a red 'X' next to each item. The widget itself shows a table with 'Weight' 120 kg 07/09/2021. The 'My Widget' section is highlighted in yellow. The 'Patient Widget' section below it shows 'No Data to Display'.

NOTE: To close overlays click outside of the overlay window