CARE4 MODULE NAME: PCS/EDM/SUR Temporary Locate a Patient

Overview

Under certain circumstances, such as when an admitted patient requires surgery or an inpatient goes to an outpatient location for treatment (i.e. Dialysis, Chemo) but will be returning to their inpatient location, the patient will need to be temporarily located. This will ensure that the patient's current location is known and lab labels print to the correct area.

Temp locating a patient is completed by the Unit receiving the patient. Example OR, Dialysis, Chemo Unit. Can be completed by Clerks/Nurses.

How to Temp Locate a Patient

Select Enter/Edit Temporary Location from the main menu

Clinical	•
Ambulatory	•
Administrative	•
Ancillary	•
Info Systems	•
Message/Task System	
RVH Reports	•
Change Your User's PIN	
Enter/Edit Temporary Location	

Identify the patient - Enter Account number or Last name, first name

Identify Patient

Select the Location from the drop-down list

	Temporary Location	
Location		-
Exp Return Date/Time		
Exp Return Mins		

Note that the patient header will now show that the patient is temporarily located by displaying (T) behind the status of ADM IN

E Return To	PCS Status Board (ROV TEST - TEST)
Care4-Rx	/h,Maria3
55 F 2	20/01/1966
ADM IN (T) V.OR

When the patient returns to their home unit it is the receiving unit that is responsible for removing the temporary location.

Follow the same steps as above to identify the patient and then remove the location.

